

Laboratory Manager Approval: Mary K. Bowman 08/19/2021
QA Manager Approval: Jeffrey Moore 08/19/2021

Glassware Maintenance

Access to this SOP shall be available within the laboratory for reference purposes; the official copy of this SOP resides on the official Georgia EPD website at <https://epd.georgia.gov/about-us/epd-laboratory-operations>. Printed copies of this SOP will contain a watermark indicating the copy is an uncontrolled copy.

1 Scope and Application

- 1.1 This is a procedure used to clean laboratory glassware. While this procedure describes a general cleanup technique, unique samples may require additional cleanup approaches to achieve interference-free glassware. *Note: Acetone should never be used in cleaning glassware to be used for TO-11A analysis. Also, TO-11A glassware should never be baked higher than 80°C. See Section 10.3 for manual cleaning of TO-11A glassware.*

2 Definitions

- 2.1 Refer to Section 3 and Section 4 of the Georgia EPD Laboratory Quality Assurance Manual for quality control definitions.

3 Interferences

- 3.1 Glassware must be scrupulously cleaned with hot water and detergent followed by de-ionized water then rinsed with methanol followed by acetone. The use of high purity reagents and solvents helps to minimize interference problems.

4 Safety

- 4.1 Refer to the Laboratory Chemical Hygiene Plan, online revision.

5 Apparatus and Equipment

- 5.1 Industrial Automatic Washer: Steris Reliance 400XLS or equivalent
5.2 Oven: Fisher Isotemp or equivalent, 105°C
5.3 Cleaning brushes: various sizes
5.4 Gloves: latex, elbow length

6 Reagents and Standards

- 6.1 Liquid Detergent
 - 6.1.1 For manual cleaning with neutral pH: Contrad, Steris Liqui-Jet 2 or equivalent
 - 6.1.2 For automated cleaning with neutral pH: Steris Glass Klenz or equivalent
- 6.2 Methanol: organic residue analysis grade
- 6.3 Acetone: organic residue analysis grade
- 6.4 De-Ionized Water or equivalent
- 7 Sample Collection**
 - 7.1 Not Applicable
- 8 Calibration**
 - 8.1 Not Applicable
- 9 Quality Control**
 - 9.1 Not Applicable
- 10 Procedure**
 - 10.1 Manual Washing (excluding TO-11A Glassware, see Sec. 10.3)
 - 10.1.1 Place small and delicate items in the sink, add a small amount of detergent (no more than 1mL) then fill with warm water and soak for 1 hour.
 - 10.1.2 Wash with a brush then rinse thoroughly with de-ionized water.
 - 10.1.3 Allow to air dry.
 - 10.1.3.1 Optional: Glassware may be placed in the glassware drying oven at 105°C and allowed to dry.
 - 10.1.3.2 Once glassware has dried in the oven, remove glassware with provided oven mitts and allow the glassware to cool to room temperature.
 - 10.1.4 Go to step 10.2.13.
 - 10.2 Automatic Washing (excluding TO-11A Glassware, see Sec. 10.3)
 - 10.2.1 Using the touchscreen on the right panel of the washer, press the “Wake” icon.
 - 10.2.2 Then press the “open door” icon to the lower left of the touchscreen to open the door. The “open door” icon looks like an “up” arrow.
 - 10.2.3 Place the glassware washing accessories to be used onto the washer manifold ensuring that washing accessory’s lower aperture aligns with the manifold port to facilitate water flow through the aperture.
 - 10.2.4 Place the appropriate racks or general purpose baskets into the washer and load the dirty glassware securely.
 - 10.2.5 Use baskets or racks if needed to cover and secure glassware and other small glassware accessories to prevent loss and breakage.
 - 10.2.6 Before starting the wash cycle, check to ensure there is adequate liquid detergent in the container beside the washer to be used in the wash cycle.
 - 10.2.6.1 If there isn’t enough liquid detergent, remove the detergent probe from the detergent container and open a new detergent container and place the detergent probe into the new detergent container.

- 10.2.7 Check to make sure that no glassware is blocking the washer door.
- 10.2.8 Press the “close door” icon at the lower left of the touchscreen to close the washer door. The “close door” icon looks like a “down” arrow.
- 10.2.9 Press the “Start Cycle” icon on the touchscreen.
- 10.2.10 Allow the washer to complete its wash cycle which includes a pre-wash, rinse, wash, rinse and drying cycles. The total wash cycle lasts approximately one hour.
- 10.2.11 After the wash cycle has completed, the touchscreen will display a “Cycle Completed” message. Press the “open door” icon and remove the cleaned glassware then press the “close door” icon.
- 10.2.12 Place the cleaned glassware onto the drying racks on the center counter top in the washroom to further air dry.
- 10.2.12.1 Optional: Glassware may be placed in the glassware drying oven at 105°C and allowed to dry.
- 10.2.12.2 Once glassware has dried in the oven, remove glassware with provided oven mitts and allow the glassware to cool to room temperature.
- 10.2.13 Once the glassware is dry and at room temperature, it is rinsed thoroughly with methanol and then rinsed again with acetone.
- 10.2.14 Allow the rinsing solvents to evaporate thoroughly in the washroom fume hood and then return the glassware to their proper storage locations.
- 10.3 Manual Washing of TO-11A Glassware
- 10.3.1 TO-11A glassware must not be rinsed with acetone and may not be baked at 105°C.
- 10.3.2 TO-11A glassware is manually washed in the same room as the TO-11A extraction fume hood.
- 10.3.3 Place the TO-11A glassware in the sink, add a small amount of detergent (no more than 1mL) then fill with warm water and soak for 1 hour.
- 10.3.4 Wash with a brush then rinse thoroughly with de-ionized water.
- 10.3.5 Allow to air dry.
- 10.3.6 Rinse glassware thoroughly with acetonitrile and allow the acetonitrile to evaporate completely in the TO-11A extraction fume hood.
- 10.3.7 Once the acetonitrile has evaporated, the TO-11A glassware may be stored in an appropriate location.

11 Calculations

- 11.1 Not Applicable

12 Waste Management

- 12.1 See GA EPD Laboratory SOP – EPD Laboratory Waste Management Standard Operating Procedures, online revision.

13 References

- 13.1 Manual for the Certification of Laboratories Analyzing Drinking Water, Criteria and Procedures Quality Assurance, Fifth Edition, Table IV-1
- 13.2 Technical Assistance Document for the National Air Toxics Trends Stations Program, Rev. 3, October 2016, Section 4.3.9.1.2

14 Reporting Limits (RLs), Precision and Accuracy Criteria and Quality Control Approach

- 14.1 Not Applicable

15 Associated LabWorks Test Codes

- 15.1 Not Applicable

Uncontrolled Copy