

Georgia Department of Natural Resources
Environmental Protection Division Laboratory

Effective Date: 06/10/2021

SOP 7-010 Rev. 4

Page 1 of 7

Laboratory Manager Approval: Ralph Schulz / 08/24/2021

QA Manager Approval: Jeffrey Moore / 08/24/2021

BALANCE USE, MAINTENANCE, AND TRAINING IN GCMS LABORATORY

Access to this SOP shall be available within the laboratory for reference purposes; the official copy of this SOP resides on the official Georgia EPD website at <https://epd.georgia.gov/about-us/epd-laboratory-operations>. Printed copies of this SOP will contain a watermark indicating the copy is an uncontrolled copy.

1. Scope and Application

- 1.1 This SOP will provide guidelines in the use of balances. These instruments are delicate and susceptible to errors in measurement if improperly calibrated, leveled, and maintained. The following document will provide the analyst with proper use of ASTM Class 1 weights, rules for validating a measurement, and preventative maintenance. Attached is a Checklist for training analysts in the use of balances (see Table 1). A copy of this Checklist should be completed and filed in the analyst's training records.

2. Definitions

- 2.1 ASTM Class 1 Weights: Weight sets that fall within the acceptable tolerance established by ASTM and NVLAP. There are 5 classes of weights: Ultra Class, Class 1, Class 2, Class 3, and Class 4. These provide a wider tolerance (and greater inaccuracy) with each higher number. Ultra Class has the smallest tolerances and is the most accurate.
- 2.2 Spirit Level: A clear usually round chamber attached to the base of a balance. This window has another circle embossed in the glass. The chamber is filled with liquid, usually ethylene glycol, and has a bubble in the liquid. The diameter of bubble will be slightly smaller than the embossed inner circle. If the balance is level, the bubble will be entirely contained within the embossed inner circle in the window of the chamber.
- 2.3 Draft Chamber: A Plexiglas™ box built around the weighing pan of the balance. It usually has several sliding doors for access to the weighing surface. These doors should be closed when a measurement is taken or a tare is made.
- 2.4 Tare: When a "tare" is made of an object, usually a beaker or other holding vessel, the balance will weigh the object and subtract from the total, giving a final value of zero on the balance display. Most tare switches are a button with a "T" "Re-zeroed" marking.
- 2.5 Zeroing the Balance: A balance should be reset, "zeroed," between each measurement. This usually means pressing the "tare" button or lever with an

empty weighing container on the pan, or the empty pan itself. This phrase is interchangeable with section 2.4, Tare.

3. **Interferences**

- 3.1 Drafts and dirt are the primary interferences for accurate weight measurements. Keep both to a minimum when weighing.
- 3.2 Weighing volatile liquids is also prone to inaccuracies due to evaporation. Speed is of the essence in this case. A set time for each measurement will give a more uniform procedure. For example, if weighing acetone into a vial, the weight will continually drop due to evaporation. If a time of say, 5 seconds, is given before the measurement is taken and the vial immediately sealed, for each measurement, the precision between multiple measurements will be good.

4. **Safety**

Refer to the Laboratory Chemical Hygiene Plan.

5. **Apparatus and Equipment**

- 5.1 Mettler PB4002 Pan Balance $\pm 0.01\text{g}$.
- 5.2 Mettler AT 200 Electronic Balance $\pm 0.0001\text{g}$
- 5.3 Set of certified ASTM Class 1 weights

6. **Reagents**

- 6.1 Not Applicable

7. **Sample Collection**

- 7.1 Not Applicable

8. **Calibration**

- 8.1.1 Calibration of a Balance will be made by a certified vendor, annually. If a balance is in immediate need of recalibration, follow the instruction in the owner's manual. If the balance can be proven, with Class 1 weights, to be accurate, it is available for use. However a re-certification of the balance should be made as soon as possible by the vendor

9. **Quality Control**

- 9.1 Use of Class 1 Weights:-Class 1 weights are expensive and easily rendered useless with improper usage. NEVER touch a Class 1 weights without gloves. The oils on your finger will cause an accumulation of dirt and change the mass of the weight, destroying its usefulness. If gloves are not used, you must use tongs provided when moving weights from the case to the balance pan. DO NOT DROP a class 1 weight. It will also change the mass of the weight and render it useless.
- 9.2 Bracketing a weight amount – If a sample, standard, or other substance has a weight taken, it must be bracketed with two Class 1 weights, one less than and one greater than the mass value. For example, if the average weight measured for sample is 50 grams and a standard is 2-5 grams, use 1 gram and a 100 gram

Class 1 weight each day to validate calibration. Do Not Use a 50-gram Class 1 weight for 50-gram weighings. Always use a higher value to bracket.

- 9.3 Balance Leveling - All Balances will have a bubble (spirit) level attached to the base. These levels will have a circle slightly larger than the bubble etched in the window. Adjust the foot on each corner of the balance until the bubble is inside the circle. Optimally, it should not touch the edge of the circle but as long as the entirety of the bubble is within the circle it is considered "level". This sounds easy but in reality it can be frustrating. The actions of a level are a mirror image of the tilt of the balance base. If you carefully raise the balance front or rear corners with your hands and tilt the base until the bubble is properly set, it gives you an idea of which foot to adjust up or down to level it.
- 9.4 Balance Cleaning – Periodically, a balance pan should be cleaned. Unplug or turn off the balance. Gently lift the pan off of the load cell. Never remove the pan if the balance is in weighing mode. Remove any dirt or debris under the pan, being careful not to let any fall in the internal parts of the balance. The top of the balance can be cleaned with a mild solvent such as methanol or propanol. Never use a harsh or abrasive cleaner on the stainless steel surface. If it is scratched to badly it will rust, and the oxide will transfer to a sample or Class 1 weight during use. Place the pan gently back of the load cell. Make sure that it does not touch anything such as a cover or part of the base. This will cause erroneous readings during use. Turn the balance back on (plug it back in) and wait 15 minutes for the balance to warm up and equilibrate.
- 9.5 Balance Placement – In order to have consistent weighings, a balance should be placed in a draft free environment. Unfortunately, this can be a problem when weighing hazardous material in a fume hood. If the balance has a draft box over the pan, the problem is less, as long as the doors are closed for each weighing. With an open pan balance, the only option is to be consistent when weighing. All balances have some type of indication that the weight is stable. It may be the appearance of a "g" next to the weight, a circle, or some other notation.
- 9.6 On a daily basis when the laboratory is open, the Mettler PB4002 uses the weights listed in Table 2 to validate the balance for that day. If a weight is out of tolerance and no problem with the balance is found then the balance is red-tagged and a service call must be placed to repair the balance.
- 9.7 On a daily basis when the laboratory is open the Mettler AT200 use the weights listed on Table 2 to validate the balance for that day. On the first day of the week the entire range of weights is used. For the remainder of the week only the 1.0g, 5.0g, 20.0g and 100g are used. If a weight is out of tolerance and no problem with the balance is found, then the balance is red-tagged and a service call must be placed to repair the balance.

10. Procedure

10.1 Mettler PB4002

- 10.1.1 The balance is turned on 24/7. If it is off, turn it on and let it warm up for 30 minutes. An overweight will display "—". Weight range is 0.1gm to 4100 gm.
- 10.1.2 To zero the balance pan press "O/T", when ready the display is "0.00" and a "g" will appear to the right of the numbers.

- 10.1.3 To weigh an object place it on the pan and wait for the 'g' to appear then record the number.
- 10.1.4 To tare a weighing container place the container on the balance pan, press "O/T" and wait for "0.00" and 'g' to appear.
- 10.1.4.1 Transfer the approximate amount of material into the weighing container and wait for the 'g' to appear then record the number.
- 10.2 Mettler AT200
- 10.2.1 The balance is turned on 24/7. If it is off, turn it on and let it warm up for 30 minutes. If 'Standby' is displayed the instrument is in the standby mode, simply press "Re-zero" and it is ready for use. An overweight will display "-----". Weight range is .01gm to 100gm. If "AUTOCALIN" appears on the display finish the weighing and step the balance to autocal in, allow balance to work for the next 15 minutes.
- 10.2.2 The balance pan is unstable if a large circle to the left of the number's readout appears. The balance is displaying a true weight value if a small flashing circle to the right of the number's readout appears. The balance has an automatic draft door.
- 10.2.1 To zero the balance pan press "Re-Zero", when it is finished it will beep, open the draft door and display "0.0000".
- 10.2.2 To weigh an object press "Re-Zero" which will open the door. Place the object on the balance and press "Print". The door will close, the balance will beep and the readout will display the weight for 7 seconds with the small circle flashing to the right of the numbers. The draft door will open automatically.
- 10.2.3 To tare a weighing container press "Re-Zero" to open the door. Place the container on the balance and press "Re-Zero" again. The balance will close the door, beep, open the door, and display "0.0000".
- 10.2.3.1 Press "Re-Zero" to open the door then transfer the amount of material into the container then press "Print" to close the door. The balance will beep, and the readout will display the weight for 7 seconds with the small flashing circle displayed to the right of the numbers. The door will open automatically.

11. Calculation

- 11.1 Not Applicable

12 References

- 12.1 Mettler Owners manual

Table 1
Balance Training Form

Analyst Name: _____

Date: _____

Person responsible for training: _____

| Question | Yes | No | Comments |
|---|-----|----|----------|
| Has the analyst read the SOP? | | | |
| Does the analyst know where the balance calibration check log is stored? | | | |
| Does the analyst know where the ASTM Class 1 weight set is stored? | | | |
| Has the analyst demonstrated the proper handling of the Class 1 weight set? | | | |
| Has the analyst demonstrated leveling a balance? | | | |
| Has the analyst demonstrated the use of balance controls? | | | |
| Has the analyst demonstrated how to clean a dirty balance? | | | |
| Has the analyst demonstrated bracketing a measurement with two Class 1 weights? | | | |
| | | | |

Analyst Signature: _____

Date: _____

Trainer's Signature: _____

Date: _____

Table 2
Balance SOP
Appendix 1
Balance Tolerance

| Mettler Model AT200 | |
|------------------------------------|------------------|
| (Weighing range 0.01-100 grams) | |
| Class S Weight, g | Tolerance |
| 0.010 | ±0.002 |
| 0.100 | ±0.002 |
| 0.200 | ±0.002 |
| 0.500 | ±0.005 |
| 1.00 | ±0.005 |
| 5.00 | ±0.005 |
| 20.0 | ±0.010 |
| 100 | ±0.010 |
| Mettler Model PB4002 | |
| (Weighing range 0.1 to 4100 grams) | |
| Class S Weight, g | Tolerance |
| 0.1 | ±0.01 |
| 1.0 | ±0.01 |
| 5.00 | ±0.01 |
| 20.0 | ±0.05 |
| 100 | ±0.1 |
| 200 | ±0.1 |
| 500 | ±0.2 |
| 1000 | ±0.3 |

**Georgia Department of Natural Resources
Environmental Protection Division**

Routine Standards Verification Form

(Class 1 Weights)

Lab. ID # _____ Balance ID# _____ Weight Set Name _____

GaDNR Balance Form -03

[illegible]

GaDNR Balance Form -03